

Product Description of Spunlace Nonwoven Polyester Viscose Nonwoven Roll For Wet Tissue


000	Spunlace Nonwoven Polyester Viscose Nonwoven Roll For Wet Tissue
000 000	Viscose + Polyester
Nonwoven Technology	Spunlace, Cross lapping
	A grade
Dotted Design	سهل
000000 00000	2.2m, could be cut into small width
:مميزات	The product has good toughness and soft feel, breathable, light, environmentally friendly, comfortable, safe and durable.
	 Medical categories: brain cotton tablets, disinfection tablets, medical dressings and so on. Sanitary categories: sanitary napkins, care mats, baby diapers, facial towels, wipes, disposable underwear, adult diapers, etc. Life categories: towels, square, compressed towels, wipes and so on. Beauty categories: cotton, cotton remover, beauty mask and so on.

Product Show of Spunlace Nonwoven Polyester Viscose Nonwoven Roll For Wet Tissue

microfiber spunlace non woven fabric custom








Related Products of Nonwoven Fabric

viscose non woven factory



Product Packaging & Shipping of Spunlace Nonwoven Polyester Viscose Nonwoven Roll For Wet Tissue



ملف الشركة

GUANGZHOU JUNQIAN NONWOVEN CO., LTD. is an enterprise located in Guangzhou city, Guangdong province,

covering an area of 58500 square meters. Concurrently, it owns its special sales company - FOSHAN GUIDE TEXTILE

CO., LTD., which is located in Foshan city. We are engaging in manufacture and sales of non woven fabric and

non woven related products, since 2005. Up to now, we have professional elite sales team, several advanced non woven

production lines with width of 1.6M, 1.9M, 2.1M, 2.5M, 3.2M, etc. Our products are well sold both in domestic and

overseas markets, widely used in areas of agriculture and gardening, medical and hygiene, personal care and beauty,

industry application, packaging and wrapping, foundation and building, hotel and restaurant, labour protection and

household, decoration, etc.

We specialized in manufacture of Spunlace Nonwoven Polyester Viscose Nonwoven Roll For Wet Tissue. نستطيع

not only control the quality and production progress, but also provide high-quality, satisfying and competitive products

for customers.

We have strict quality management system, and equipped with professional non woven test. Our products meet various

of standards and own various of certificates and reports, such as SGS, MSDS, BUREAU VERITAS, CE, Nemko,

TüV Rheinland, ATC, CTI, etc. We can design and produce according to customer's requirements and specifications,

depending on the advanced technology and excellent management. Relying on many clients' trust, we are also offering

cousultant and agency services on non woven products. With the business principle of "Quality First, Reputation Builds

Future", we keep on diligently and excelsior working, to improve marketing competitive abilities, and finally accomplish

company goal.



FAQ of Spunlace Nonwoven Polyester Viscose Nonwoven Roll For Wet Tissue

0000	
Q1: Do you have it in stock ?	A1: Yes,we have.
Q2: What is your MOQ ?	A2: Customized:1ton
Q3: How can I get the sample ?	A3: If you need sample to test, we can make it as per your request. But we need to charge sample fees. If it is our regular product in stock, you just pay freight cost and sample is free.
Q4: What is your package detail ?	A4: individually packed for each pcs,or customer's requirements
Q5: What should I do if I cannot find what I looking for on your website ?	A5: Please feel free to contact us by email and send detailed information for the product you look for.We will check whether it is available.

Q6: What about lead time for mass production ?	A6: Honestly, it depends on the order quantity and the season you place order. Generally speaking,the lead time is about 2-4 weeks. So we suggest that you start inquiry as earlier as possible.
Q7: Is it all right to make customer's own brand name ?	A7: 000000 00000 0000 .00000000 OEM. 0 000000 000000 000000 000 0 000000000
000 00000 00000 000 :8 0 000000	A8: 000 .00000000 000000 000 000 000 000 0