

## Product Description of Two-component Non Woven Fabric For Medical And Restaurant Applications JQht7050-w2-90

Item	<a href="#">Two-component Polyester Non Woven Fabric</a> For Medical And Restaurant Applications JQht7050-w2-90
Raw Material	Two-component
Nonwoven Technic	Spunbond
Grade	A grade
Design	Plain
Colors	White Colour
Standard Width	2.2m, could be cut into small width
Weight	as request are available
Feature	Non-toxic, non-irritating, product stability, high strength
Applications	It widely use in surgical coats, <a href="#">face masks</a> ; sanitary napkins for women, baby diapers, adult masts, protective suits and overalls, agricultural bumper, geotextiles, cables and motor insulators, battery separators, <a href="#">packaging materials</a> , garment linings, cloth, waterproof material and so on.

## Product Show of Two-component Non Woven Fabric For Medical And Restaurant Applications JQht7050-w2-90

A close-up photograph of a nonwoven fabric material. The fabric has a textured, fibrous appearance with a scalloped or wavy edge on the left side. The color is a light, neutral tone. Two website URLs are overlaid on the image in a semi-transparent font.

[www.nonwovenproductsulipper.com](http://www.nonwovenproductsulipper.com)

[www.guldenonwoven.com](http://www.guldenonwoven.com)

A close-up photograph showing several layers of nonwoven fabric. The top layer is a smooth, light-colored material. Below it, there are several layers of a more fibrous, textured material, likely a filter or separator. The layers are slightly offset, creating a stepped effect. The background is a dense, fibrous material, possibly the same as the layers shown.

[www.nonwovenproductsulipper.com](http://www.nonwovenproductsulipper.com)

[www.guidenonwoven.com](http://www.guidenonwoven.com)



Product Packaging and Shipping of Two-component Non Woven Fabric For  
Medical And Restaurant Applications JQht7050-w2-90



Packed in rolls,  
with 2 " or 3 "  
core inside,  
strong poly  
bags outside



transporting the goods  
to the port



loaded in container

## FAQ of Two-component Non Woven Fabric For Medical And Restaurant Applications JQht7050-w2-90

Question	Answer
Q1: Do you have it in stock ?	A1: Yes,we have.
Q2: What is your MOQ ?	A2: customized:5 ton
Q3: How can I get the sample ?	A3: If you need sample to test, we can make it as per your request. But we need to charge sample fees. If it is our regular product in stock, you just pay freight cost and sample is free.
Q4: What is your package detail ?	A4: individually packed for each pcs,or customer's requirements

Q5: What should I do if I cannot find what I looking for on your website ?	A5: Please feel free to contact us by email and send detailed information for the product you look for.We will check whether it is available.
Q6: What about lead time for mass production ?	A6: Honestly, it depends on the order quantity and the season you place order. Generally speaking,the lead time is about 3-6 weeks. So we suggest that you start inquiry as earlier as possible.
Q7: Is it all right to make customer's own brand name ?	A7: Sure. We warmly welcome the OEM clients. We can discuss the package, design, etc.
Q8: When can I get the price ?	A8: Usually we quote within 24 hours after we get your inquiry detail. If you are very urgent to get the price, please call us or tell us in your email so that we will quote your the price as soon as possible.